## APM 440 plus

The APM 440 plus alternating pressure system is suited for the treatment of pressure ulcers up to grade IV (according to EPUAP).

The control unit is equipped with a digital sensor control, which regulates the inner pressure of the mattress, creating a perfect surface pressure. Again, the use of two compressors guarantees safety in case of a compressor malfunction. Adjustment is accomplished with the aid of a touch pad surface, equipped with visually easy to comprehend pictograms. The result is a quick and easy operation.

- I bi-elastic, breathable PU/PES tubes, connection coupling wi safety lock, external CPR valve
- I sensor operated pressure mar between cycles switch off, au key lock, alternating pressure, and care mode, 3 time cycles 5/10/15 minute selection
- I holds up to 230 kg patient weight

















Accessories: integrable pillow

# **Mattress Specifications**

Item no. 991500X

ontrol unit

Measurements: approx. 39 x 21 x 12 cm Weight: approx. 5 kg

2 membrane compressors, digital control with pressure sensor, alternating pressure, static mode, care mode, 3 time cycles (5/10/15 min.), key lock function, low pressure alarm,

power failure alarm, back-up function

Mattress

Measurements inflated incl.

base mattress: approx. 200 x 90 x 22 cm Weight: approx. 12.5 kg

20 cells, 3 head cells static, non-bending connecting tubes\*, self-closing connection couplings, CPR emergency ventilation

#### Materials

Cells: nylon with polyurethane coating Cover: bi-elastic PU/PES, breathable,

slip-retardant
Base: nylon cover with integrated

### Included in delivery

transport bag for mattress and control unit, foam mattress base with nylon cover

pocket for 2 cm base mattress

## Pressure ulcer grade:

• I • II • III • I

Max. patient weight: 230 kg

\* incl. transportation safety lock