



ACTIVE WHEELCHAIRS

AKTIV / IRIS

KURY spol. s r.o.
Zdechovice 28, 533 11, Czech Republic
tel.: +420 466 936 139 • fax: +420 466 936 138
e-mail: kury@kury.cz • www.kury.cz

KURY

IRIS X1

AKTIV X1

Innovated version of the popular active ultralight wheelchair with **rigid frame**:

- dynamic and attractive design for highly active users
- tough compact aluminum frame with excellent driving characteristics
- individual combinations and adjustment possibilities
- ergonomics and folding backrest
- adjustment of the angle of the backrest and footrest
- straps for multipoint adjustment of backrest stiffness
- adjustment of center of gravity
- total weight from 10.1 kg

IRIS X1 design model **NATURAL***


- fresh design with clear varnish frame and footrest
- rear wheels Dynamic DYS
- carbon mudguards
- front wheels Al-PU 100 x 32
- grey textile backrest and seat combined with fabric DOBBY®

Ultralight **folding** wheelchair for highly active users:

- attractive design, high performance
- durable folding aluminum and tough frame
- 12 positions of adjustment of center of gravity
- backrest with straps for stiffness adjust
- seat tension adjustment
- textile backrest combined with fabric DOBBY®
- positioning footrest
- total weight from 11.6 kg

AKTIV X1 - trend*

- ANTRACIT color of frame and cross fork and rear wheel holder - colored with ELOX according to the offer
- rear wheels ENERGY
- ABS or carbon mudguards
- front wheels Al-PU (100) 125 x 32
- black textile backrest and seat combined with fabric DOBBY®



attractive, high-quality fabric cover combined with fabric DOBBY®



adjustment of center of gravity



straps for multipoint adjustment of backrest stiffness



adjustment of the angle of the backrest and its complete folding



height-adjustable and angle adjustable footrest



attractive, high-quality fabric cover combined with fabric DOBBY®



seat tension adjustment



straps for multipoint adjustment of backrest stiffness




12 positions of adjustment of center of gravity




positioning and folding footrest

* Product includes optional equipment incl. additional cost items.



carbon classic



white carbon



ABS plastic

mudguards



black



silver



blue



red


colored forks



AKTIV S

Ultralight mechanical folding wheelchair for active users:

- attractive design
- durable folding aluminum and tough frame
- easy folding and transportability
- 12 positions of adjustment of center of gravity
- positioning and folding footrest
- a wide range of settings and accessories
- total weight from 11.8 kg



adjustment of center of gravity



positioning and folding footrest



seat tension adjustment



solid construction and functional design



completely folding



IRIS



KURY

All the wheelchairs are made individually to measure of user. You may find all the input parameters, dimensions and complete accessories incl. order codes in the manufacturer's questionnaire (see also website www.kury.cz), or please contact KURY company directly. The manufacturer reserves the right to modify technical and other parameters without prior notice. KURY company's products are certified and produced in Europe.
Made in Czech Republic. © KURY spol. s r.o. 11/2013