

boris. As clean as can be.

GMFCS
Level III-V

kg
max. 150



26-46 cm

31-48 cm

(-15° - 30°)



I am your comfortable shower and toilet chair. Whether cleaning your teeth or showering, my seat tilt can be adapted to any situation. My anatomically shaped seat surface gives you the required support, stabilises your upper body and can be adjusted in height. You can sit even more comfortably on my insertion pads. These can easily be removed and cleaned. Because I have been produced from highly robust material, you and your parents will have a lot of fun with me for a long time to come. Did you know that I am available in five sizes and can carry up to 150 kg?



You can sit really comfortably on the **removable, optional seat pad.**



Swing a way arm supports also help to stabilise you



**Seat pan**

Anatomically moulded, provides stability in the pelvic area and stabilises the upper body.

Headrest

can be adjusted in height and angle.

Table

with padded armrest.

Insertion pad

prevents pressure marks on the body.

Robust sub-frame

made from stainless steel.

Footplate

can simply be folded up for easier access.


Foot pedal

adjusts the seat pan continuously from (-)5° to 30°.

**Toilet equipment**

You can use me independently of a toilet as a toilet chair.

I am also available with a larger hygiene opening in sizes 2-4.





Data	Size 0	Size 1	Size 2	Size 3	Size 4
Seat width (rear)	26 cm	31 cm	36 cm	41 cm	46 cm
Seat depth	31 cm	36 cm	41 cm	46 cm	48 cm
Seat shell width (front)	31 cm	36 cm	41 cm	46 cm	51 cm
Back height	49 cm	51 cm	56 cm	61 cm	65 cm
Seat angle	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°
Can be positioned over a toilet	47 - 57 cm	45 - 55 cm	42 - 52 cm	42 - 52 cm	42 - 52 cm
Total width	56 cm	61 cm	63 cm	71,5 cm	74 cm
Total length	65 cm	65 cm	65 cm	70 cm	70 cm
Max. lower leg length	38 cm			48 cm	
Front edge of hygiene opening to rear edge of sub-frame	47 cm	49 cm	51 cm	54 cm	57 cm
Front opening width	A	10 cm	10 cm	10 cm	10 cm
Rear opening width	b	12 cm	12 cm	13,5 cm	13,5 cm
Rear opening width*	c	-	-	18,5 cm	18,5 cm
max. load	60 kg	80 kg	100 kg	100 kg	150 kg
Weight	28 kg	30 kg	33 kg	35 kg	38 kg
HMV No.	33.40.04.1014				

* Only applies to sizes 2 b, 3 b and 4 b



You can find more information on accessories & harnesses in the price list.



Scan this QR code and watch videos about the product directly on our homepage.

